PrO-Xtract

By Centurion Scientific Ltd









CE marked ISO 9001 regulated and now IVD certified too

Welcome to our new PrO-Xtract range.

Designed with the busy Chef in mind

Derived from our Long running LED and LCD display centrifuge spanning over 25 years of Centrifuge manufacturing.

Simplicity of use, is paramount and places this emphasis accordingly. All are supplied with a standard rotor or a combination to suit your needs.

Ease of Use

Pre programmed with a variety of menus or choose your own, then just press start to go, once finished, a buzzer will sound and you press Lid open to retrive your samples, simplicity itself. We offer both Ambient and Refrigerated version, for your convienience.

Larger models are available from our huge range if required.


Please call us

See our User, Standard and Safety Features on the page opposite.

Built to last, strong construction yet offering a sleek modern design that will fit into any modern laboratory.


No compromises on quality, only the best components are used in the manufacturing process.

By working with all our suppliers, our reliability has been retained and improved to give all our range's a **3 year warranty** as standard.



Centurion Scientific Ltd

PrO-Xtract Chef Centrifuge



Built to last, strong construction yet offering a sleek modern design that will fit into any modern laboratory.

No compromise on quality only the best components are used in the manufacturing process.

By working with all our suppliers, our reliability has been retained and improved to give all our range a 3 year warranty as standard.

User features to PrO-Xtract	Advantages
Large LCD display with touch screen control	Clarity & ease of use
Rotor Recognition	Safe selection of rotors
10 Accelerations rates 10 deceleration rates	Sample assistance
108 program memory (12 pages of 9)	Multi department & Users
Timer 0-9999 minutes & Hold in 1 secon increments	Total flexibility
Pulse Short run	Fast Pelleting
Run in Speed (rpm) or rcf (G) in 1 rpm increments	Accuracy
W2t Plots tim, accel rpm and decel	True sample repeatability
Sound <60db (rotor dependent)	Quiet
Orientation acceleration rate	Prevents initial sample side deposits

Standard to PrO-Xtract	Advantages
Extra thick Stainless steel bowl	Easy Clean & rust free
Port to Lid	Tachometer speed certification
Alloy & Steel frame with composite's	Strong Light & Quiet
World Leading Industrial grade inverter	Reliability & Strength
Best Quality European Brushless motor	Quiet, cool & reliable
High technology airflow (ambient mode)	Cooler running
All Centrifuge have user accesible service sections to access all saftey parameters	

Saftey features to PrO-Xtract	Advantages
Multi point lid locking	Lid Safety
Emergency lid release	Power failure release
Lid gas struts	Lid dropping safety
Lid lock detection	To run Lid must be shut
Imbalance detection	Easy "eye" volume loading
Overspeed Sensor	Safe detection of speed
Set inverter values	Electronic safety of speed
Barrier ring	Extra metal protection of chamber
Motor overheat sensor	Safe motor protection

Refrigerated PrO-Xtract	Advantages
Variable -9°C to +40°C with stand by cooling	Flexibility (see rotors)
CFC free R404A gas	Ozone friendly
Valve control of Gas	Efficient temperature control
PID control of valve	Supreme temperature accuracy

Centurion Scientific Ltd Centrifuges comply to all relevant EU standards of Quality and medical Devices IEC 61010 and CE conformity test marks Emission, immunity to EN/IEC 61326-1, Class B

Index

- Page 5 PrO-Xtract 7. Home Chef
 Ambient Centrifuge 400ml capacity
 4000 Rpm, 2,200 G
- Page 6 PrO-Xtract5 and PrO-Xtract 5R
 Ambient or Refrigerated 1 Litre Centrifuge
 6000 Rpm, 4,800G
- **Page 7 -** Rotors available to Pro-Xtract 5/5R 4 x 250ml Swing out rotor 6 x 100ml fixed angle rotor
- Page 8 PrO-Xtract 3 and PrO-Xtract 3R
 Ambient or Refrigerated 3 Litre Centrifuge
 10000 Rpm, 15,650 G
- **Page 9 -** Rotors available to Pro-Xtract 3/3R 4 x 750ml Swing out rotor 6 x 250ml fixed angle rotor
- Page 10 Endorsements

 Meet the Chefs already exploring the possibilities
- Page 11 Accessories

 Available Bottles & Consumables

PrO-Xtract7. Home Chef **400ml Capacity**



Display indicative only

Speed	500-4000 Rpm
	(1 Rpm steps)
Rcf max	2,200 G
Timer	0-9999 Mins & Hold
	(1 sec steps)
DimensionsHWD	275 x 325 x 470mm
Weight	21.8 Kg (without rotor)
Power	120 watts
Memory	108 programs
Accel rates	10
Decl rates	10

Swing Out Rotor (Included)



BRK5505

Rotor type Tube size max Minimum Speed Maximum Speed Maximum Rcf (G) Radius max cms Sample tube angle (°) Acceleration time (secs) Deceleration time (secs) Autoclaveable (frequency

Rotors

4 x 100ml 45 x 100mm 500 Rpm 4,000 Rpm 2200 12.5 0 30 30 121°C (20)

Bottles available see page 11

PrO-Xtract5 Ambient Centrifuge 1 Litre Capacity



Display indicative only

Speed	500-6000 Rpm
	(1 Rpm steps)
Rcf max	4,800G
Timer	0-9999 Mins & Hold
	(1 sec steps)
DimensionsHWD	310 x 400 x 500mm
Weight	31.6 Kg (without rotor)
Power	309 watts
Memory	108 programs
Accel rates	10
Decl rates	10


PrO-Xtract5R Refrigerated Centrifuge 1 Litre Capacity



Display indicative only

Speed	500-6000 Rpm
	(1 Rpm steps)
Rcf max	4,800G
Timer	0-9999 Mins & Hold
	(1 sec steps)
DimensionsHWD	315 x 450 x 635mm
Weight	62.2 Kg (without rotor)
Power	448 watts
Memory	108 programs
Accel rates	10
Decl rates	10
Temp	9°C to +40°C PID
	controlled to +/- 1°C

PrO-Xtract 5/R Available Rotors



121°C (20)

Bottles available - see page 11

121°C (20)

Rotors

Rotor type

Tube size max Minimum Speed

Maximum Speed

Maximum Rcf (G)

Radius max cms

Sample tube angle (°)

Acceleration time (secs)

Deceleration time (secs)

Autoclaveable (frequency

Buckets B5250 (Set 4)

Sealed Lids B5419

PrO-Xtract3 Ambient Centrifuge 3 Litre Capacity



Display indicative only

Speed	500-10000 Rpm
	(1 Rpm steps)
Rcf max	22,000 G
Timer	0-9999 Mins & H old
	(1 sec steps)
DimensionsHWD	375 x 600 x 630mm
Weight	63.5 Kg (without rotor)
Power	750 watts
Memory	108 programs
Accel rates	10
Decl rates	10

PrO-Xtract3R Refrigerated Centrifuge 3 Litre Capacity



Display indicative only

Speed	500-10000 Rpm
	(1 Rpm steps)
Rcf max	22,000 G
Timer	0-9999 Mins & Hold
	(1 sec steps)
DimensionsHWD	410 x 998 x 630mm
Weight	120 Kg (without rotor)
Power	1200 watts
Memory	108 programs
Accel rates	10
Decl rates	10
Temp	9°C to +40°C PID
	controlled to +/- 1°C

PrO-Xtract 3/R Available Rotors



Bottles available - see page 11

Rotors

Rotor type

Tube size max

Minimum Speed

Maximum Speed

Maximum Rcf (G)

Sample tube angle (°)

Acceleration time (secs)

Deceleration time (secs)

Autoclaveable (frequency

Radius max cms

Buckets B3000 (Set 4)

Sealed Lids B5319

Endorsements

The Pro Xtract range has been well received by its numerous users. The simplicity and flexibility of the range has ensured that Chefs from many different training backgrounds and establishments have been able to introduce recipes inspired by this unique piece of equipment. By breaking down ingredients in to its singular component's Chefs are able to create very distinct flavours by using ingredients differently.

Aulis - Simon Rogan

Adjacent to the 2 Michelin Star restaurant L'Enclume based in Cartmel, Aulis offers guests an exclusive experience of personally cooked dishes by the development Chefs. Used as a development kitchen, the team investigate new techniques and technologies for use in Simon's many restaurants. Our equipment is used on a regular basis and Simon has enjoyed finding new ways to use ingredients and exploring the various layers to a simple ingredient.

Claridges-Simon Rogan

Further to the success of the Centrifuge in Aulis when Simon took over Claridges in the early part of 2014 it was ensured that a larger PrO-Xtract 3 would be an integral part of the Kitchen. Chefs ensure that recipes that are developed in Aulis are produced on a much larger scale for diners of Fera.

On Air Dining - Jamie Cairns

Jamie was very excited to receive his Centrifuge after meeting with us at Hotelympia. On Air Dining provides an exclusive catering service for Private Charter flights. Based in London On Air Dining services flights from Gatwick, Luton, Stansted and many private airfields. Guests of this service are now experiencing dishes inspired by using our centrifuges like clear Consommé and Fruit oils. All the dishes are prepared in his Kitchen and supplied with step by step instructions to recreate these beautiful dishes at 35,000 feet.

South Lodge Hotel - Matt Gillan

Set in beautiful surroundings South Lodge Hotel in West Sussex offers guests Michelin Star dining in 'The Pass'. Matt knew in principle the function of a centrifuge but was unaware of how easy and readily available they were to the catering industry. After demonstrating the Pro-Xtract range to Matt he was sold instantly on how this could add to the dining experience. Using the Centrifuge for making Broad Bean Butter is Matt's favourite dish, he uses it to add a twist to simple starters.

Bottles available 100-750ml

Safety

Please be aware, not any plastic container is suitable, as the G force involved will split or rupture the vessel and could be damaging to the centrifuge.

Molecular Menu

Although we will kepp a list of updated menu's to try on our Global website.

If you would like to contribute, we will publish in confidence or in detail.

(Free advertising for you)

Further Information

We trust that the information contained in this Sales literature proves informative.

However if you have any questions, or require a larger capacity, please do not hesitate tocontact us. Centurion is a manufacturer of a huge range of centrifuges and can accomodate most of our customers needs



To review all our ranges please visit our global website

www.centurionscientificglobal.com

